TMX[®] PARA RÜCKENTRIGGER

TMX[®] PARA

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PHYSIO@TMX-TRIGGER.DE

TMX[®] PARA BACK TRIGGER

TMX[®] PARA OUR MOST EFFECTIVE SPINAL TRIGGER

THE INTELLIGENT SOLUTION:

With its horizontal double trigger heads the **TMX[®] PARA** is able to trigger your pain areas to the left and right of your spine simultaneously.

WHY TRIGGER THE LEFT AND RIGHT OF THE SPINE AT THE SAME TIME?

Starting from our rear pelvis, humans have an ascending musculature which extends to the right and left of our spine. Deeper muscle layers serve the inner stability of our vertebral bodies. The superficial muscle layers serve more to the mobility of our individual spinal column sections.

Severe and long-term pain is very common in these areas.

This persistent, often „deep“ back pain frequently extends over a larger area of the back. The pain is difficult to localize and is oftentimes referred to as „diffuse deep“ back pain.

WHY IS TMX[®] PARA CALLED PARA?

The word PARA is derived from the Latin name *paravertebral*, which means „left and right parallel to the spine“.

THE TMX[®] 3-STEPS METHOD: APPLICATION

1. Place the **TMX[®] PARA** on the ground. Locate your pain point and lie down on the **PARA**. **IMPORTANT:** Make sure that your spine is between the two trigger heads! Now put pressure on the painful area and try to relax despite the pain. Remain in this position for at least 30 seconds.

2. Now start mobilizing by moving your arms or legs up and down. Meanwhile, inhale and exhale calmly several times. Alternate between mobilizing and the resting starting position every 30 seconds. Do this switch for 3-4 minutes.

3. Finally, stretch your back.

3 STEPS TO A FREE MUSCLE

1. TRIGGER

2. MOBILIZE

3. STRETCH

TMX[®] PARA – 3 DIFFERENT TRIGGER VARIANTS:

1. **ORIGINAL**
 - head diameter 3 cm
 - superficial muscle groups (low intensity)
2. **PLUS**
 - head diameter 2 cm
 - deep muscle groups (medium intensity)
3. **TWIN**
 - head diameter 2 x 0,8 cm
 - deep muscle groups (high intensity)

TRIGGERS IN ACUTE PAIN:


- A** pain area; trigger right and left of the spine;
- B** surrounding area which ideally should also be triggered afterwards

TRIGGER FOR DIFFUSE PAIN OR AS PREVENTION:

- C** In case the pain area cannot be localized, or simply as means of regeneration, trigger along the spine from top to bottom.

TRIGGERING PROMOTES REGENERATION:

- D** Trigger the area of your thoracic spine to regulate the stress center (sympathetic nervous system). Trigger your posterior cranial area and the entire sacrum to activate the relaxation nerve (parasympathetic nervous system).

AREAS OF APPLICATION:

- pain in the lumbar spine
- pain between lower ribs and iliac crest
- pain in the thoracic spine
- pain between the shoulder blades
- general back pain

NOTES

IMPORTANT: Depending on the tension and „gluing“ of the myofascial tissue, a multi-day application of up to a week may be necessary. In case of sore muscles on the following day, a rest day should be observed between the applications. To prevent recurring tension, the application should be continued even if there are no pain symptoms.

DO NOT USE WITH

Signs of inflammation (warmth, swelling, redness), after operations, acute herniated disc, fractures, metal in the spine, osteoporosis, swirl gliding (for further information visit our website)